

At a glance

B. Braun hygiene products for instruments and surfaces

Product	Area of application													Scope of efficacy							Active disinfectant ingredient										Application								
	Surface disinfection	Surface cleaning	Inventory	Large surfaces	Small surfaces/fast disinfection	Sanitation of suction systems and amalgam separators	Instrument disinfection	Instrument cleaning	Disinfection of hemodialysis machines	Bedpan preparation	For metal, glass, ceramic	For thermo-sensitive items	For flexible endoscopes	Bacteria	Tuberculosis bacteria	MRSA	Spores	Fungi	Limitedly virucidal (incl. HBV, HCV, HIV)	Limitedly virucidal PLUS	Virucidal	HBV/HIV	Noroviruses (NOV)	Alcohols	Aldehydes	Formaldehyde-free	Phenol derivatives	Quarternary ammonium compounds	Alkylamines	Peracetic acid	Citric acid	Active chlorine	Ready-to-use	pH of ready-to-use solution	Concentrate to be diluted	Concentration for application (DGHM/VAH)	Exposure time (DGHM/VAH)	Concentration for application / exposure time (IFSG Section 18)	
Helipur®							●	●			●			●	●	●		●	●			●				●									9.5	●	1.5% 3.0%	1 h 5 mins	6% 2-6 h
Helipur® H plus N							●				●	●	●	●	●	●		●	●		●	●	●		●	●									5.0	●	1.0% 1.5%	30 mins 15 mins	
Stabimed® Ultra							●	●			●	●	●	●	●	●	●		●	●		●	●			●									8.0	●	2.0%	15 mins	
Stabimed® fresh							●	●			●	●	●	●	●	●		●	●				●					●							9.0	●	0.5% 1.0%	15 mins 5 mins	
Helizyme								●			●	●	●													●									7.0	●	1.0%	5 mins	
Cleaner N								●			●	●	●													●									7.0	●	1.0%	15 mins	
Helimatic® Rinse neutral NF								●			●															●									7.0	●	0.05- 0.2%		
Helimatic® Cleaner neutral								●			●	●														●									7.0	●	0.1- 0.5%		
Helimatic® Cleaner alkaline								●		●	●	●														●									11.5	●	0.3- 0.8%		
Helimatic® Neutralizer C								●			●															●									3.0	●	0.05- 0.3%		
Helimatic® Cleaner MA								●			●	●	●													●									10.0	●	0.2- 1.0%		
Helimatic® Cleaner enzymatic NF								●			●	●	●													●									7.0	●	0.1- 0.5%		
Helimatic® Disinfectant							●				●	●	●	●	●	●		●	●		●	●	●		●										7.0	●	1.0%	5 mins	
Helimatic® Latriniser								●		●	●															●									5.0	●	0.3- 1.0%		
Tiutol® dent for dentistry						●								●	●	●		●	●		●	●	●			●								12.5	●	3.0%	60 mins		
Meliseptol® HBV Wipes	●		●		●									●	●	●		●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	7.0		Conc.	1 mins		
Meliseptol® Wipes sensitive	●	●	●		●									●	●	●		●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	7.0		Conc.	1 mins		
Meliseptol® rapid	●		●		●									●	●	●		●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	7.0		Conc.	1 mins		
Meliseptol® Foam pure	●	●	●		●									●	●	●		●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	7.0		Conc.	1 mins		
Meliseptol®	●		●		●									●	●	●		●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	7.0		Conc.	1 mins		
Hexaquant® XL	●	●	●	●	●									●	●	●		●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	9.0	●	2.0% 1.0%	5 mins 15 mins			
Hexaquant® pure	●	●	●	●	●									●	●	●		●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	8.5	●	0.5% 1.0%	15 mins 5 mins			
Melsept® SF	●	●	●	●										●	●	●		●	●	●	●	●	●	●	●	●	●	●	●	●	●	5.0	●	0.5% 2.0%	1 h 15 mins				
Meliseptol® Wipes ultra	●	●	●		●									●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	11.0		Conc.	2 mins			
Citric acid solution 50%								●						●	●			●	●		●	●	●			●							2.0	●	3.0%	15 mins 83°C			
Tiutol® KF								●						●	●			●	●		●	●	●			●							12.5	●	3.0%	15 mins 60°C			