



B. Braun disinfectants and cleaning agents for manual or automated preparation

Manual instrument cleaning and disinfection at a glance	Thermostable materials	Thermo-sensitive materials	Flexible endoscopes	Cleaning	Effective against: * Bacteria	* Fungi	* Spores	Limitedly virucidal (incl. HBV, HCV, HIV)	Virucidal
Helipur®	✓			✓	✓	✓		✓	
Helipur® H plus N	✓	✓	✓		✓	✓		✓	✓
Stabimed® Ultra	✓	✓	✓	✓	✓	✓	✓	✓	✓
Stabimed® fresh	✓	✓		✓	✓	✓		✓	
Helizyme	✓	✓	✓	✓					

Automated instrument cleaning and disinfection at a glance	Instrument disinfection	Instrument cleaning	Bedpan preparation	Rinsing	Neutralizing
Helimatic® Cleaner alkaline		✓			
Helimatic® Cleaner MA		✓			
Helimatic® Cleaner enzymatic NF		✓			
Helimatic® Cleaner neutral		✓			
Helimatic® Neutralizer C					✓
Helimatic® Disinfectant	✓				
Helimatic® Rinse neutral NF				✓	
Helimatic® Latriniser			✓	✓	