

Compatibility of surface and instrument disinfection

		Surfaces						Instruments						
		Melsept® SF	Hexaquart® XL	Meliseptol® HBV Wipes	Meliseptol® acute	Meliseptol® Foam V	Meliseptol® Wipes sensitive/ acute	Hexaquart® pure	Helipur® H plus N	Stabimed® fresh	Helipur®	Stabimed® Ultra	Cleaner N	Helizyme®
Surfaces	Melsept® SF	✓	✗	✓	✓	✓	✓	✓	✓	✗	✓	✓	✓	✓
	Hexaquart® XL	✗	✓	✓	✓	✓	✓	✓	✗	✓	✓	✓	✓	✓
	Meliseptol® HBV Wipes	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Meliseptol® acute	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Meliseptol® Foam V	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Meliseptol® Wipes sensitive/ acute	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Hexaquart® pure	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Instruments	Helipur® H plus N	✓	✗	✓	✓	✓	✓	✓	✓	✗	✓	✓	✓	✓
	Stabimed® fresh	✗	✓	✓	✓	✓	✓	✓	✗	✓	✓	✓	✓	✓
	Helipur®	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Stabimed® Ultra	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Cleaner N	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Helizyme®	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓